

Conjugate Release Pad

Document No.: MSDS/1002

Revision Date: 02 January 2022

1. Product & Company Identification

Product Name Conjugate Release Pad

Product Information Conjugate Release Pad is critical to Lateral Flow Assay

> performance, the Synthetic Pad must dry the Conjugate without damage and release the conjugate rapidly when the analyte

flow to the Pad.

These Pad can be used for all kind of lateral flow testes like monitoring, ovulation, detecting infectious disease organisms, analysing drugs of abuse and measuring other analyst important to human physiology hCG, Cancer, Whole Blood test (HIV, Malaria) and also agricultural applications etc. using

sandwich, competitive or inhibition assays.

Product Code

Cat. No.	Thickness (mm)	Material Type	Weight (mg/cm²)	Wicking Time (Sec.)
16613	0.42	Untreated Poly.	100	
16613T01	0.42	Treated Poly.	100	30
16615	0.53	Untreated Poly.	135	
16615T01	0.53	Treated Poly.	135	30
16613N	0.40	Untreated Poly.	100	
16613NT01	0.40	Treated Poly.	100	25

Company Name : **Axiflow Biotech Private Limited**

Factory Address:- Plot No. 328, sector56, Phase-V, EHTP, HSIIDC, Kundli, Sonipat, Haryana, 131028 INDIA **Address** :

Phone No. +91-11-40578224 / 25

E-mail arvind@axiflow.co.in, jatin@axiflow.co.in :

Website www.axiflow.co.in

2. Physical Properties

Physical State Solid

Colour White

Physical Form Sheet/Reel

Odour Natural Polyester Odour

Boiling Point Not applicable

3. Chemical Properties

Water Solubility : Insoluble

pH : Neutral

Melting Point : Not applicableEvaporation Rate : Not applicable

Volatility : Not applicable

4. Hazard Identification

Major Health Hazards : No significant target effects reported.

Physical Hazards: Flammable solid. May cause eye & skin irritation.

Precautionary Measures : Do not eat, drink or smoke when using this product. Avoid

contact with eyes, skin and clothing. Keep away from heat,

sparks and flame. Keep container closed. Wash thoroughly

after handling.

Potential Health Effect : Ingestion - Although, Unlikely to be ingested. It may cause

irritation of digestive system and kidney damage.

Skin Contact - No known any allergic reaction reported.

Eye Contact - Although, Unlikely to happened, It may cause

irritation and burns in eyes.

5. First Aid Measures

Eye Contact : It is unlikely that emergency treatment will be required.

Wash with large amounts of water or normal saline until no

evidence of chemical remains. Get medical attention

immediately if needed.

Skin Contact: No skin problem has been reported. If it occurs, wash skin with

soap or mild detergent & large amounts of water,

Seek medical attention immediately if needed.

Inhalation : No Inhalation hazard. Still if any irritation arises, move

exposed person to fresh air. If breathing is irregular or if

respiratory arrest occurs, provide artificial respiration or

oxygen by trained personnel, Get medical help immediately.

Ingestion: Wash out mouth with water. Do not induce vomiting unless

directed to do so by medical personnel. Never give anything by

mouth to an unconscious person. Get medical attention

immediately.

6. Fire Fighting Measures

Flammability : Flammable solid. Runoff to sewer may create fire or explosion

hazard.

Extinguishing media : Dry chemical, sand, CO₂, regular foam, water

Specific Hazards : Promptly isolate the scene by removing all persons from the

vicinity of the incident if there is a fire. No action shall be taken

involving any personal risk or without suitable. training

Move containers from fire area if this can be done without risk

Use water spray to keep fire exposed containers cool.

Special protective equipment for fire-fighters

Fire-fighters should wear appropriate equipment and self-contained breathing apparatus (SCBA) with a full face piece operated in positive pressure mode.

7. Accidental Release Measures

Methods for Clean up : Move containers from spill area. Vacuum or sweep up material

and place in a designated, labeled waste container. Use spark

proof tools or explosion-proof equipment.

Personal Precautions : No action shall be taken involving any personal risk or without

suitable training. Evacuate surrounding areas. Do not touch or

walk through spilled material. Shut off ignition sources. No

flares, smoking or flames in hazard area. Provide adequate

ventilation. Wear appropriate respirator when ventilation is

inadequate.

Environmental Precautions: Avoid dispersal of spilled material and run-off and contact

with soil, waterways, drains and sewers.

8. Handling and Storage

Handling : Keep the product in original packing. Avoid direct sunlight and

organic solvents. Keep away from heat, sparks, flame and

excessive humidity or moisture. Eating, drinking and smoking

should be prohibited in areas where this material is handled,

stored and processed. Use explosion proof electrical equipment.

Storage : Store in a dry and well-ventilated area place. Store below 30 °C.

It is very good if the temperature is below 20 °C. Keep container

slightly closed to prevent moisture absorption and

contamination.

9. Exposure Control/Personal Protection

Ventilation : Based on available information, additional ventilation is not

required. General ventilation is adequate.

Personal Protection Equipment

Eye Protection : Eye protection not required under normal conditions

Face Protection: Face protection not required under normal conditions

Skin Protection: Protective gloves not required under normal conditions

Respiratory Protection: No respirator is required under normal conditions

10. Stability & Reactivity

Reactivity : Stable at normal temperatures and pressure.

Stability : Stable at normal temperatures and pressure.

Condition to avoid : Avoid heat, flames, sparks and other sources of ignition.

Incompatibilities : Acids, alcohols, amine, bases, ketones, oxidizing materials

Hazardous: Thermal decomposition products or contact with acids:

Decomposition : Oxides of nitrogen and oxides of carbon.

11. Toxicological Information

Routes of Entry : Inhalation, Ingestion

Primary Effects : Harmful if inhaled or ingested.

On the skin : Severely irritates skin and may result in severe pain.

On the eye : Causes eye irritation.

Ingestion : Harmful if swallowed. Swallowing may result in severe

stomach pain.

Inhalation : May be harmful if inhaled. May cause irritation of the

respiratory tract and mucous membranes

12. Ecological Information

Aquatic Toxicity : Toxicity is expected to be low based on insolubility in water.

13. Disposal Considerations

Waste Disposal

The generation of waste should be avoided or minimized wherever possible. Significant quantities of waste product residues should not be disposed of via the foul sewer but processed in a suitable effluent treatment plant. Dispose of surplus and non-recyclable product via a licensed waste disposable contractor.

Incineration and landfill should only be considered when recycling material and it is not feasible. This container must be disposed of in a safe way. Cares should be taken when handling emptied containers that have not been cleaned or rinsed out. Empty containers may retain some product residues. Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers.

14. Transport Information

Transport Information

These products are classified as "Flammable Solid" in the

meaning of DOT Classification.

15. Regulatory Information

HCS Classification

Flammable solid, target organs effects

Protective Equipment : Gloves, Lab Coat, Dust respirator. Wear appropriate respirator

when ventilation is inadequate. Splash goggles.

16. Other Information

This Material Safety Data Sheet serves to complete but not to replace the technical product data sheets. The information and recommendations set forth above are taken from sources believed to be accurate as of the date hereof; however, Axiflow Biotech makes no warranty with respect to the accuracy of the information or the suitability of the recommendations and assumes no liability to any user thereof. The information contained in this sheet does not constitute a hazard assessment and should not be used in place of the user's own assessment of workplace risks as required by other health and safety legislation. The user bears sole liability for the precautions required when using the product. The regulatory texts indicated herein are intended to aid the user to fulfill his obligations. This list is not to be considered complete and exhaustive. It is the user's responsibility to ensure that he is subject to no other obligations than those mentioned. The above information is believed to be correct, but does not purport to be all-inclusive and shall be used only as a guide.

Axiflow Biotech shall not be held liable for any damage resulting from handling or from contact with the above product.

MANUFACTURER DISCLAIMER

Information given herein is offered in good faith as accurate and experience of supplier but without guarantee. Conditions of use and suitability of the product for particular uses are beyond our control; all risks of use of the product are therefore assumed by the user. Nothing is intended as a recommendation for uses.